



# Recombinant Escherichia coli Electron transport complex protein RnfC (rnfC)

<b>Product Code</b>	CSB-MP542656ENX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1LEQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain SMS-3-5 / SECEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLKLFSAFRK NKIWDFNGGI HPPEMKTQSN GTPLRQVPLA QRFVIPLKQH IGAEGELCVS VGDKVLRGQP LTRGRGKMLP VHAPTSQTVT AIAPHSTAHP SALAEHSVII DADGEDCWIP RDGWDDYRSR SREELIERIH QFGVAGLGGA GFPTGVKLQG GGDKIETLII NAAECEPYIT ADDRMLQDCA AQVVEGIRIL AHILQPREIL IGIEDNKPQA ISMLRAVLAD SHDISLRVIP TKYPSGGAKQ LTYILTQKQV PHGGRSSDIG VLMQNVGTAY AVKRAVIDGE PITERVVTLT GEAIARPGNV WARLGPVVRH LLNDAGFCPS ADQMVMGGP LMGFTLPWLD VPVVKITNCL LAPSANELGE PQEEQNCIRC SACADACPAD LLPQQLYWFS KGQQHDKATT HNIADCIECG ACAWVCPSNI PLVQYFRQEK AEIAAIRQEE KRAAEAKARF EARQARLERE KAARLERHKS AAVQPAKDK DAIAAALARV KEKQAQATQP IVIKAGERPD NSAIIAAREA RKAQARAKQA ELQQTNDAAT VADPRKTAVE AAIARAKARK LEQQQANAEP EQQVDPRKAA VEAAIARAKA RKLEQQQANA EPEQQVDPRK AAVEAAIARA KARKLEQQQA NAEPEEQVDP RKA AVAAAIA RVQAKKAAQK KVVNED
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsxC
<b>Protein Names</b>	Recommended name: Electron transport complex protein RnfC
<b>Expression Region</b>	1-676
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.