



Recombinant Escherichia coli UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase (arnB)

Product Code	CSB-BP542978ENT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1X8W6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / DH10B)
Purity	>85% (SDS-PAGE)
Sequence	MSEFLPFSRP AMGVEELAAV KEVLESGWIT TGPKNQALEQ AFCQLTGNQH AIAVSSATAG MHITLMALKI GKGDEVITPS LTWVSTLNMI SLLGATPVMV DVDRDTLMVT PEAIESAITP RTKAIIPVHY AGAPADIDAI RAIGERYGIA VIEDAAHAVG TYYKGRHIGA KGTAIFSFA IKNITCAEGG LIVTDNENLA RQLRMLKFHG LGVDAYDRQT WGRAPQAEVL TPGYKYNLTD INAAIALTQL VKLEHLNTRR REIAQQYQQA LAALPFQPLS LPAWPHVHAW HLFIIIRVDEQ RCGISRDALM EALKERGIGT GLHFRAAHTQ KYRERFPTL SLPNTEWNSE RICSLPLFPD MTTADADHVI TALQQLAGQ
Source	Baculovirus
Target Names	arnB
Protein Names	Recommended name: UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase EC= 2.6.1.87 Alternative name(s): UDP-(beta-L-threo-pentapyranosyl-4'-ulose diphosphate) aminotransferase Short name= UDP-Ara4O aminotransferase UDP
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.