



# Recombinant Escherichia coli Pyrimidine-specific ribonucleoside hydrolase rihB (rihB)

<b>Product Code</b>	CSB-MP543648ENT
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1X860
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12 / DH10B)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEKRKIILDC DPGHDDAIAI MMAAKHPAID LLGITIVAGN QTLDKTLING LNVCKLEIN VPHYAGMPQP IMRQQIVADN IHGETGLDGP VFEPLTRQAE STHAVKYIID TLMASDGDIT LVPVGPLSNI AVAMRMQPAI LPKIREIVLM GGAYGTGNFT PSAEFNIFAD PEAARVVFTS GVPLVMMGLD LTNQTVCTPD VIARMERAGG PAGELFSDIM NFKTKTQFEN YGLAGGPVHD ATCIGYLINP DGIKTQEMYV EVDVNSGPCY GRTVCDELGV LGKPANTKVG ITIDTDWFWG LVEECVIRGYI KTH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rihB
<b>Protein Names</b>	Recommended name: Pyrimidine-specific ribonucleoside hydrolase rihB EC=3.2.2.8 Alternative name(s): Cytidine/uridine-specific hydrolase
<b>Expression Region</b>	1-313
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.