



Recombinant Escherichia coli D-cysteine desulfhydrase (dcyD)

Product Code	CSB-MP544476ENX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1LQT7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SMS-3-5 / SECEC)
Purity	>85% (SDS-PAGE)
Sequence	MPLHNLTRFP RLEFIGAFTP LEYLPRFSDY LGREIFIKRD DVTPMAMGGN KLRKLEFLAA DALREGADTL ITAGAIQSNH VRQTAAVAAK LGLHCVALLE NPIGTTAENY LTNGNRLLLD LFNTQIEMCD ALDTPNAQLE ELATRVEAQQ FRPYVIPVGG SNALGALGYV ESALEIAQQC EGAVNISSVV VASGSAGTHA GLAVGLEHLM PESELIGVTV SRSVADQLPK VVNLQQAIAK EELTASAEI LLWDDYFAPG YGVPNDEGME AVKLLARLEG ILLDPVYTGK AMAGLIDGIS QKRFKDEGPI LFIHTGGAPA LFAYHHPHV
Source	Mammalian cell
Target Names	dcyD
Protein Names	Recommended name: D-cysteine desulfhydrase EC= 4.4.1.15
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.