



# Recombinant Escherichia coli GMP reductase (guaC)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP544846ENX                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                       |
| <b>Uniprot No.</b>       | B1LG40                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen Species</b> | Escherichia coli (strain SMS-3-5 / SECEC)                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Sequence</b>          | <p> MRIEEDLKLKLG FKDVLRPKR STLKSRSDVE LERQFTFKHS GQSWSGVPPII<br/> AANMDTVGTF SMASALASFD ILTAVHKHYS VEEWQAFINN SSADVLKHVM<br/> VSTGTSDADF EKTKQILDNL PALNFVCIDV ANGYSEHFVQ FVAKAREAWP<br/> TKTICAGNVV TGEMCEELIL SGADIVKVGI GPGSVCTTRV KTGVGYPQLS<br/> AVIECADAHAH GLGGMIVSDG GCTTPGDVAK AFGGGADFVM LGGMLAGHEE<br/> SGGRIVEENG EKFMLFYGMS SESAMKRHVG CVAEYRAAEG KTVKLPLRGP<br/> VENTARDILG GLRSACTYVG ASRLKELTKR TTFIRVQEQQE NRIFNNL </p> |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Target Names</b>      | guaC                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Names</b>     | Recommended name: GMP reductase EC= 1.7.1.7 Alternative name(s): Guanosine 5'-monophosphate oxidoreductase Short name= Guanosine monophosphate reductase                                                                                                                                                                                                                                                                                 |
| <b>Expression Region</b> | 1-347                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                      |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                              |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                       |