



Recombinant Danio rerio Nucleoredoxin (nxn)

Product Code	CSB-EP016224DIL
Abbreviation	nxn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q503L9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSEFLVLLG ERLVNGEKAE VDVQALGSRL SLLGLYFGCS LNGPCKQFNA SLTEFYSKFK KSSEHKDKLE IVFISSDQDQ KQWQDFLQEM QWPALPFKDR HKMKLWNKY KVTSLPSLVF IDAATGKVVC RGNLLVVRDD PKGLEFPWGP KPF AEV VSGP LLRNNRQTTD STALEGSYVG VYFSAHWCPP CRSLTRVLVE SYRKVKETGQ KFEIVFVSAD RSEESFTQYF SEMPWLAVPY SDEARRSRLN RLYGIQGIPT LILLDTEGHM ITRQGRVEIL NDPDCGLFPW HPRPVLELSE SNAVQLHEGP CLVLFVDAEE EGELDPAKEL IQPIAEKIMA KYKAKEEETP LLFFVAGEDD MSDSLRDYTN LPEAAPLLTI LDMSARAKYV KDVEEITPAV VEQFVSGFLA EKLKPEPI
Source	E.coli
Target Names	nxn
Protein Names	Recommended name: Nucleoredoxin EC= 1.8.1.8
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.