



Recombinant Bubalus bubalis Leptin (LEP)

Product Code	CSB-MP012870BWX
Abbreviation	LEP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5J732
Product Type	Recombinant Protein
Immunogen Species	Bubalus bubalis (Domestic water buffalo)
Purity	≥85% (SDS-PAGE)
Sequence	VPIRKVQDD TKTLIKTIVT RINDISHTQS VSSKQRVTGL DFIPGLHPLL SLSKMDQTLA IYQQILTSLP SRNVVQISND LENLRDLLHL LAASKSCPLP QVRALESLES LGVVLEASLY STEVVALSRL QGSLQDMLRQ LDLSPGC
Source	Mammalian cell
Target Names	LEP
Protein Names	Recommended name: Leptin Alternative name(s): Obesity factor
Expression Region	22-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.