



# Recombinant Human Midnolin (MIDN)

<b>Product Code</b>	CSB-MP013822HU
<b>Abbreviation</b>	MIDN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q504T8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPQPGGARS CRRGAPGGAC ELGPAAEAAP MSLAIHSTTG TRYDLAVPPD ETVEGLRKRL SQRLKVPKER LALLHKDTRL SSGKLQEFVGDGSKLTLVP TVEAGLMSQA SRPEQSVMQA LESLTETQVS DFLSGRSPLT LALRVGDHMM FVQLQLAAQH APLQHRHVLA AAAAAAARG DPSIASPVSS PCRPVSSAAR VPPVPTSPSP ASPSPITAGS FRSHAASTTC PEQMDCSPTA SSSASPGAST TSTPGASPAP RSRKPGAVIE SFVNHAPGVF SGTFSGLTLP NCQDSSGRPR RDIGTILQIL NDLLSATRHY QGMPPSLAQL RCHAQCSPAS PAPDLAPRTT SCEKLTAAAPS ASLLQGQSQI RMCKPPGDRL RQTENRATRC KVERLQLLLQ QKRLRRKARR DARGPYHWSP SRKAGRSDSS SSGGGGSPSE ASGLGLDFED SVWKPEVNPDIKSEFVVA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MIDN
<b>Protein Names</b>	Recommended name: Midnolin Alternative name(s): Midbrain nucleolar protein
<b>Expression Region</b>	1-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.