



Recombinant Human Interferon-induced protein 44-like (IFI44L)

Product Code	CSB-EP011015HU-B
Abbreviation	IFI44L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53G44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVTTRLTWN DENHLRKLLG NVLSLLYKS SVHGGSIEDM VERCSRQGCT ITMAYIDYNM IVAFMLGNYI NLHESSTEPN DSLWFSLQKK NDTTEIETLL LNTAPKIIDE QLCVRLSKTD IFIICRDNKI YLDKMITRNL KLRFYGHRQY LECEVFRVEG IKDNLDDIKR IIKAREHRNR LLADIRDYRP YADLVSEIRI LLVGPVGS GK SFFNSVKSI FHGHVTGQAV VGSDITSITE RYRIYSVKDG KNGKSLPFML CDTMGLDGAE GAGLCMDDIP HILKGCM PDR YQFN SRKPIT PEHSTFITSP SLKDRIHCVA YVLDINSIDN LYSKMLAKVK QVHKEVLNCG IAYVALLTKV DDCSEVLQDN FLNMSRSMTS QSRVMNVHKM LGIPISNILM VGNYASDLEL DPMKDILILS ALRQMLRAAD DFLEDLPLEE TGAIERALQP CI
Source	E.coli
Target Names	IFI44L
Protein Names	Recommended name: Interferon-induced protein 44-like
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.