



Recombinant Rat Smad nuclear interacting protein 1 (Snip1)

Product Code	CSB-MP021931RA
Abbreviation	Snip1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5M9G6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MKAGKSERER SSRRRHRSRGD ALATVVVKQE RLSPEPVAHR RPDAPAASPP PPAAESGSAG HRGSRARGAS RSPAKKKSKS SGRRSKSPRT KRSRSPHYST VKVKQEREDH PRRGREDRQH RELSEQEHRR ARNSERDRHR GHARQRRSSD ERPVSGQGRD RDSQILQAQE EERDFNNARR REHRQQNESA GAEAQEVIPR PAGNKNKEVP VKEKPSFELS GALLEDTNTF RGVVIKYSEP PEARIPKKRW RLYPFKNDEV LPVMIYHRQS AYLLGRHRR ADIPIDHPSC SKQHAFVQYR LVEYTRADGT VGRRVKPYII DLGSGNGTFL NNKRIEPQRY YELKEKDVLK FGFSSREYVL LHESSDTSEL DRKEDEDEEE EEEMVSDSS
Source	Mammalian cell
Target Names	Snip1
Protein Names	Recommended name: Smad nuclear interacting protein 1
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.