



Recombinant Pelodiscus sinensis Glyceraldehyde-3-phosphate dehydrogenase (GAPDH)

Product Code	CSB-EP009232PEL-B
Abbreviation	GAPDH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R2J2
Product Type	Recombinant Protein
Immunogen Species	Pelodiscus sinensis (Chinese softshell turtle) (Trionyx sinensis)
Purity	>85% (SDS-PAGE)
Sequence	MVKIGVNGFG RIGRLVTRAA FTCDKVQIVA INDPFIDLNY MVYMFKYDST HGRFHGTVKA ENGKLVINGQ AITIFQERDP ANIKWGDAGA EYVVESTGVF TTTEKAS AHL KGGAKRVVIS APSADAPMFV MGVNHEKYDN SLKVVSNASC TTNCLAPLAK VIHDNYGMVE GLMTTVHAIT ATQKTVDGPS GKLWRDGRGA AQNIIPASTG AAKAVGKVIP ELNGKLTGMA FRVPTPNVSV VDLTCRLQKP AKYDDIKKVM KAASEGPLKG ILGYTEDQVV SSDFNGDSRS SIFDAAAGIA LNDNFVKLVS WYDNEFGYSN RVVDLLVHMA SKE
Source	E.coli
Target Names	GAPDH
Protein Names	Recommended name: Glyceraldehyde-3-phosphate dehydrogenase Short name= GAPDH EC= 1.2.1.12 Alternative name(s): Peptidyl-cysteine S-nitrosylase GAPDH EC= 2.6.99.-
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.