



Recombinant Mouse Growth factor receptor-bound protein 10 (Grb10)

Product Code	CSB-EP009887MO-B
Abbreviation	Grb10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60760
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNNDINSSVE SLNSACNMQS DTDAPLLED GQHASNQGAA SSSRGQPQAS PRQKMQRSQP VHILRRLQEE DQQLRTASLP AIPNPFPELT GAAPGSPPSV APSSLPPPPS QPPAKHCGRC EKWIPGENTR GNGKRKIWRW QFPPGFQLSK LTRPGLWTKT TARFSKKQPK NQCPTDTVNP VARMPQSME KLRLRKDVKV FSEDGTSKVV EILTDMTARD LCQLLVYKSH CVDDNSWTLV EHPQLGLER CLEDEIVVQ VESTMPSESK FLFRKNIYAKY EFFKNPVNFF PDQMVNWCQQ SNGGQAQLLQ NFLNTSSCPE IQGFLQVKEV GRKSWKKLYV CLRRSGLYYS TKGTSKEPRH LQLLADLEES SIFYLIAGKK QYNAPNEHGM CIKPNKAKTE MKELRLLCAE DEQIRTCWMT AFRLLYGML LYQNYRIPQR KGLPPPFNAP MRSVSENSLV AMDFSGQIGR VIDNPAAEQS AALEEGHAWR KRSTRMNILS SQSPLHPSTL NAVIHRTQHW FHGRISREES HRIKQQGLV DGLFLLRDSQ SNPKAFVLT CHHQKIKNFQ ILPCEDDGQT FFTLDDGNTK FSDLIQLVDF YQLNKGVLPC KLKHH CIRVA L
Source	E.coli
Target Names	Grb10
Protein Names	Recommended name: Growth factor receptor-bound protein 10 Alternative name(s): GRB10 adapter protein Maternally expressed gene 1 protein
Expression Region	1-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. Overexpression of some isoforms of the encoded protein inhibits tyrosine kinase



activity and results in growth suppression. This gene is imprinted in a highly isoform- and tissue-specific manner. Alternatively spliced transcript variants encoding different isoforms have been identified.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.