



# Recombinant Human ATP synthase mitochondrial F1 complex assembly factor 1 (ATPAF1)

<b>Product Code</b>	CSB-EP002424HU-B
<b>Abbreviation</b>	ATPAF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5TC12
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GGA DSSGVGAEAE LQANPFYDRY RDKIQLLRRS DPAAFESRLE KRSEFRKQPV GHSRQGDFIK CVEQKTDALG KQSVNRGFTK DKTLLSSIFNI EMVKEKTAEE IKQIWQQYFA AKDTVYAVIP AEKFDLIWNR AQSCPTFLCA LPRREGYEFF VGQWTGTELH FTALINIQTR GEAAASQLIL YHYPELKEEK GIVLMTAEMD STFLNVAEAQ CIANQVQLFY ATDRKETYGL VETFNLRPNE FKYMSVIAEL EQSGLGAELEK CAQNQNKT
<b>Source</b>	E.coli
<b>Target Names</b>	ATPAF1
<b>Protein Names</b>	Recommended name: ATP synthase mitochondrial F1 complex assembly factor 1 Alternative name(s): ATP11 homolog
<b>Expression Region</b>	58-328
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes an assembly factor for the F(1) component of the mitochondrial ATP synthase. This protein binds specifically to the F1 beta subunit and is thought to prevent this subunit from forming nonproductive homooligomers during enzyme assembly. Alternatively spliced transcript variants have been identified, but the biological validity of some of these variants has not been determined.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.