



Recombinant Mouse TNF receptor-associated factor 3 (Traf3)

| | |
|--------------------------|---|
| Product Code | CSB-BP024148MO |
| Abbreviation | Traf3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q60803 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MESSKKMDAA GTLQPNPPLK LQPDRGAGSV LVPEQGGYKE KFKVKTVEDKY KCEKCRLVLC NPKQTECGHR FCESCMAALL SSSSPKCTAC QESIHKDKVF KDNCKREIL ALQVYCRNEG RGCAEQLTLG HLLVHLKNEC QFEELPCLRA DCKEKVLKRD LRDHVEKACK YREATCSHCK SQVPMIKLQK HEDTDCPCVV VSCPHKCSVQ TLLRSELSAH LSECVNAPST CSFKRYGCVF QGTNQQIKAH EASSAVQHVN LLKEWSNSLE KKVSLQNES VEKNKSIQSL HNQICSFEIE IERQKEMLRN NESKILHLQR VIDSQAEKLELDKEIRPFR QNWEEADSMK SSVESLQNRV TELESVDKSA GQAARNTGLL ESQLSRHDQM LSVHDIRLAD MDLRFQVLET ASYNGVLIWK IRDYKRRKQE AVMGKTLISLY SQPFYTYGYFG YKMCARVYLN GDGMGKGTHL SLFFVIMRGE YDALLPWPFFK QKVTLMMLMDQ GSSRRHLGDA FKPDPNSSSF KKPTGEMNIA SGCPVFVAQT VLENGTYIKD DTIFIKVIVD TSDLPDP |
| Source | Baculovirus |
| Target Names | Traf3 |
| Protein Names | Recommended name: TNF receptor-associated factor 3 Alternative name(s): CD40 receptor-associated factor 1 Short name= CRAF1 TRAFAMN |
| Expression Region | 1-567 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Target Details | This protein is a member of the TNF receptor associated factor (TRAF) protein family. TRAF proteins associate with, and mediate the signal transduction from, members of the TNF receptor (TNFR) superfamily. This protein participates in the signal transduction of CD40, a TNFR family member important for the activation of the immune response. This protein is found to be a critical component of the lymphotoxin-beta receptor (LTbetaR) signaling complex, |



which induces NF-kappaB activation and cell death initiated by LTbeta ligation. Epstein-Barr virus encoded latent infection membrane protein-1 (LMP1) can interact with this and several other members of the TRAF family, which may be essential for the oncogenic effects of LMP1. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.