



# Recombinant Mouse Mitogen-activated protein kinase 3 (Mapk3)

<b>Product Code</b>	CSB-EP013456MO-B
<b>Abbreviation</b>	Mapk3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q63844
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AAAAAAPGG GGGEPRGTAG VVPVVPGEVE VVKGQPFVDG PRYTQLQYIG EGAYGMVSSA YDHVRKTRVA IKKISPFEHQ TYCQRTLREI QILLRFRHEN VIGIRDILRA PTLEAMRDVY IVQDLMETDL YKLLKSQQLS NDHICYFLYQ ILRGLKYIHS ANVLHRDLKP SNLLINTTCD LKICDFGLAR IADPEHDHTG FLTEYVATRW YRAPEIMLNS KGYTKSIDIW SVGCILAEML SNRPIFPCHK YLDQLNHILG ILGSPSQEDL NCIINMKARN YLQSLPSKTK VAWAKLFPKS DSKALDLLDR MLTFNPNKRI TVEEALAHYPY LEQYYDPTDE PVAEEPFTFD MELDDLPER LKELIFQETA RFQPGAPEGP
<b>Source</b>	E.coli
<b>Target Names</b>	Mapk3
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase 3 Short name= MAP kinase 3 Short name= MAPK 3 EC= 2.7.11.24 Alternative name(s): ERT2 Extracellular signal-regulated kinase 1 Short name= ERK-1 Insulin-stimulat
<b>Expression Region</b>	2-380
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the MAP kinase family. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act in a signaling cascade that regulates various cellular processes such as proliferation, differentiation, and cell cycle progression in response to a variety of extracellular signals. This kinase is activated by upstream kinases, resulting in its translocation to the nucleus where it phosphorylates nuclear targets. Alternatively spliced transcript variants encoding different protein isoforms have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.