



Recombinant *Xenopus tropicalis* Growth factor receptor-bound protein 2 (grb2)

| | |
|--------------------------|---|
| Product Code | CSB-EP009889XBF-B |
| Abbreviation | grb2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q66I13 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEAVAKYDFK ATADDELSFK RGDILKVLNE ECDQNWYKAE LNGKDGFIK NYIEMKPHPW FFGKIPRAKA EEMLGKQRHD GAFLIRESES APGDFSLSVK FGNDVQHFKV LRDGAGKYFL WVKFNSLNE LVDYHRSTSV SRNQQIFLRD IEQVPQVHGG DRATNLLQQP TYVQALDFD PQEDGELGFR RGDFIQVVDN SDPNWWKGTC LGQTGMFPRN YVTPVNRNM |
| Source | E.coli |
| Target Names | grb2 |
| Protein Names | Recommended name: Growth factor receptor-bound protein 2 Alternative name(s): Adapter protein GRB2 SH2/SH3 adapter GRB2 |
| Expression Region | 1-229 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |