



Recombinant Pongo abelii Antizyme inhibitor 1 (AZIN1)

Product Code	CSB-BP002484PYX
Abbreviation	AZIN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R7K3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	MKGFIDDVDY SVGLLDEGTN LGNVIDNYVY EHTLTGKNAF FVGDLGKIVK KHSQWQNVVA QIKPFYTVKC NSAPAVLEIL AALGTGFACS SKNEMALVQE LGVPPENIIY ISPCQVSQI KYAAKGVNI MTCNEIELK KIARNHPNAK VLLHIATEDN IGEEGNMKF GTTLKNCRHL LECAKELDVQ IIGVKFHVSS ACKESQVYVH ALSDARCVFD MAGEIGFTMN MLDIGGGFTG TEFQLEEVNH VISPLLDVYF PEGSGVKIIS EPGSYVSSA FTLAVNIIAK KVVENDKFPS GVEKTGSDEP AFMYYMNDGV YGSFASKLSE DLNTIPEVHK KYKEDEPLFT SSLWGPSCDE LDQIVESCLL PELNVGDWLI FDNMGADSFH EPSAFNDFQR PAIYYMMSFS DWYEMQDAGI TSDSMMKNFF FVPSCIQLSQ EDSFSAEA
Source	Baculovirus
Target Names	AZIN1
Protein Names	Recommended name: Antizyme inhibitor 1 Short name= AZI Alternative name(s): Ornithine decarboxylase antizyme inhibitor
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.