



# Recombinant Pongo abelii Antizyme inhibitor 1 (AZIN1)

<b>Product Code</b>	CSB-EP002484PYX
<b>Abbreviation</b>	AZIN1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R7K3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKGFIDDVDY SVGLLDEGTN LGNVIDNYVY EHTLTGKNAF FVGDLGKIVK KHSQWQNVVA QIKPFYTVKC NSAPAVLEIL AALGTGFACS SKNEMALVQE LGVPPENIIY ISPCQVSQI KYAAKGVNI MTCDNEIELK KIARNHPNAK VLLHIATEDN IGEEGNMKF GTTLKNCRHL LECAKELDVQ IIGVKFHVSS ACKESQVYVH ALSDARCVFD MAGEIGFTMN MLDIGGGFTG TEFQLEEVNH VISPLLDVYF PEGSGVKIIS EPGSYVSSA FTLAVNIIAK KVVENDKFPS GVEKTGSDEP AFMYYMNDGV YGSFASKLSE DLNTIPEVHK KYKEDEPLFT SSLWGPSCDE LDQIVESCLL PELNVGDWLI FDNMGADSFH EPSAFNDFQR PAIYYMMSFS DWYEMQDAGI TSDSMMKNFF FVPSCIQLSQ EDSFSAEA
<b>Source</b>	E.coli
<b>Target Names</b>	AZIN1
<b>Protein Names</b>	Recommended name: Antizyme inhibitor 1 Short name= AZI Alternative name(s): Ornithine decarboxylase antizyme inhibitor
<b>Expression Region</b>	1-448
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.