



Recombinant Human PRAME family member 5 (PRAMEF5)

Product Code	CSB-EP018622HU-B
Abbreviation	PRAMEF5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TYX0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSIRTPPRL ELAGRSLLRD QALAMSTLEE LPTELPPLF MEAFSRRRCE ALKLMVQAWP FRRLPLRPLI KMPCLEAFQA VLDGLDALLT QGVHPRRWKL QVLDLQDVCE NFWMVWSEAM AHGCFLNAKR NKKPVQDCPR MRGQQPLTVF VELWLKNRTL DEYLTCLLLW VKQRKDLLHL CCKKLKILGM PFRNIRSILK MVNLDIQEV EVNCKWVLP I LTQFTPYLGH MRNLQKLVLS HMDVSRVYVSP EQKKEIVTQF TTQFLKLCCL QKLSMNSVSF LEGHLDQLLS CLKTSLKVL TITNCVLES D LKHL SQCP SI SQLKTL DLSG IRLTNYSLVP LQILLEKVA TLEYLDLDDC GIIDSQVNAI LPALSRFEL NTF SFCGNPI SMATLENLLS HTIILKNLCV ELYPAPRESY DADGTLCWSR FAQIRAELMK RVRDLRHPKR ILFCTDCCPD CGNRSFYDLE ADQCCC
Source	E.coli
Target Names	PRAMEF5
Protein Names	Recommended name: PRAME family member 5
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.