



Recombinant Human PRAME family member 8 (PRAMEF8)

Product Code	CSB-EP018625HU-B
Abbreviation	PRAMEF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VWM4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSIRAPPRL ELARQRLLRD QALAISTMEE LPRELFPTLF MEAFSRRRCE TLKTMVQAWP FTRLPLGSLM KSPHLESLSKS VLEGVDVLLT QEVRPRQSKL QVLDLRNVDE NFCDIFSGAT ASFPEALSQK QTADNCPGTG RQQPFMVFID LCLKNRTLDE CLTHLLEWVK QKRGLLHVCC QELQVFGMPI HSIIEVLNMV ELDCIQEVEV CCPWELSTLV KFAPYLGQMR NLRKLVLFNI RASACIPPDN KGQFIARFTS QFLKLDYFQN LSMHSVSFLE GHLDQLLRCL QASLEMVVMV DCLLSESDLK HSWCPSIRQ LKELDLRGVT LTHFSPEPLT GLLEQVVATL QTLDEDCGI MDSQLSAILP VLSRCSQLST FSFCGNLISM AALENLLRHT VGLSKLSLEL YPAPLESYDT QGALCWGRFA ELGAELMKTL RDLRQPKIIV FCTVPCPRCG IRASYDLEPS HCLC
Source	E.coli
Target Names	PRAMEF8
Protein Names	Recommended name: PRAME family member 8
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.