



Recombinant Chicken WAS protein family homolog 1 (WASH1)

Product Code	CSB-EP025972CH-B
Abbreviation	WASH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZKA6
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MTTVAQKHFL EGQTYSVPLI QPDLRREEAV QQVADALQYL QKVSGDIFNR ISQRVETSRA QLQAISERV LAQAKIEKIK GSKKAIKVFS SAKYPAPERL QEYCSIFAGA EDPSKQKWPR HKIQSKHRML DEKSLQEKLK YFPVCVNTKI HQEDDAEEGL GSLPRNISSL SLLLLFNTTE NLYKKYVFLD PLAGAVTKTH VALETETEEK LFDAPLSITE RGQLDRQVAE NYFYVPDLGQ VPEIDVPSYL PDLPGIAADL MYSADLPGI APSAPSSAIP ELPTFSTESV EPSQADLQDP GLLPPPPPPP PPPPPVMPTT VPPPPPLPQP TAPSEPARTA SEDSSKTVPA ASVQGAPKEV VNPSTGRASL LESIRQAGGI GKANLRVKE RKLEKKKQKE QEQVRATGQG GDLMSDLFNK LVLRRKGISG KGP GASANPD APGSPAGAF RMSDSIPPLP PPQQPPGEED EDDWES
Source	E.coli
Target Names	WASHC1
Protein Names	Recommended name: WAS protein family homolog 1
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.