



Recombinant Mouse Ral guanine nucleotide dissociation stimulator-like 2 (Rgl2)

Product Code	CSB-BP019625MO
Abbreviation	Rgl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61193
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MLPRPLRLLL DTTPPGGVVL SSFRSRDPEE GGDPGGRAVG GGQEEEEDEEE EEASVSVWDE EEDGATFTVT SRQYRPLDPL APLPPPRSSR RLRAGTLEAL VRHLLDARTA GADMMFTPAL LATHRAFTST PALFGLVADR LEALESHPPG ELERTTGVAI SVLSTWLASH PEDFGSEVKG QLDRLESFLL RTGYAAREGV VGGSADLIRN LRARVDPRAP DLPKPLALPG DSPADPTDVL VFLADHLAEQ LTL LDAELFL NLIPSQCLGG LWGHRDRPGH SHLCPSVRAT VTQFNKVAGA VVSSVLGATS IGEGPREVTV RPLRPPQRAR LLEKWIRVAE ECRLLRNFSS VYAVVSALQS SPIHRLRAAW GETTRDSL RV FSSLCQIFSE EDNYSQSREL LMQEVKQPP VEPHKKAPR SGFRGGGVVP YLGTFLKDLV MLDAASKDEL ENGYINFDKR RKEFAILSEL LRLQKECRGY DLRPNSDIQQ WLQGLQPLTE AQSHRVSCEV EPPGTSDDSPA ARTPRPTLVI TQWTEVLGSV GGPTPLVSWD RPSVGGDEVP GTPAPLLTRL AQHMKWPSVS SLDSALESSP SLHSPADPGH LSPPASSPRP SRGHRRSASC GSPLSGNTGE GTSRSAGCGG GVSGPGSSDC RIIRVQMELG EDGSVYKSIL VTSQDKAPSV ISRVLKKNNR DSAVASEFEL VQLLPGDREL TIPHSANVFY AMDGASHDFL LRQRRRPSAA TPGSHSGPSA SGTPPSEGGG GSFPRIKATG RKIARALF
Source	Baculovirus
Target Names	Rgl2
Protein Names	Recommended name: Ral guanine nucleotide dissociation stimulator-like 2 Short name= RalGDS-like 2 Alternative name(s): RalGDS-like factor Ras-associated protein RAB2L
Expression Region	1-778
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.