



# Recombinant Human Fc receptor-like B (FCRLB)

<b>Product Code</b>	CSB-BP008561HU
<b>Abbreviation</b>	FCRLB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6BAA4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATL EKPILSLHPP WTTIFKGERV TLKCDGYHPL LLELQPISTL WYLGHLLLPS HKKSIEVQTP GVYRCQTRGA PVSDPIHLSV SNDWLILQVP YAPVFEGEPL VLRCRGWYDK VVYKLHYYHD GQAVRYFHSS ANYTVLQARA SDSGRYQCSG TMRIPVESAP MFSAKVAVTV QELFRAPVLR VMGPREARGA ALGGVVLRC DTRLHPQKRDT PLQFAFYKYS RAVRRFDWGA EYTVPEPEVE ELESYWCEAA TATRSVRKRS PWLQLPGPGS PLDPASTTAP APWAAALAPG NRPLSFRKPP VSRSVPLVTS VRNTTSTGLQ FPASGAPTAG PPACAPPTPL EQSAGALKPD VDLLLREMQL LKGLLSRVVL ELKEPQALRE LRGTPTPTS HFAVSPGTPE TTPVES
<b>Source</b>	Baculovirus
<b>Target Names</b>	FCRLB
<b>Protein Names</b>	Recommended name: Fc receptor-like B Alternative name(s): Fc receptor homolog expressed in B-cells protein 2 Short name= FREB-2 Fc receptor-like and mucin-like protein 2 Fc receptor-like protein 2 Fc receptor-related protein Y
<b>Expression Region</b>	18-426
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.