



Recombinant Human Fc receptor-like B (FCRLB)

Product Code	CSB-EP008561HU
Abbreviation	FCRLB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BAA4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ATL EKPILSLHPP WTTIFKGERV TLKCDGYHPL LLELQPISTL WYLGHLLLPS HKKSIEVQTP GVYRCQTRGA PVSDPIHLSV SNDWLILQVP YAPVFEGEPL VLRRCRGWYDK VVYKLHYYHD GQAVRYFHSS ANYTVLQARA SDSGRYQCSG TMRIPVESAP MFSAKVAVTV QELFRAPVLR VMGPREARGA ALGGVVLRC DTRLHPQKRDT PLQFAFYKYS RAVRRFDWGA EYTVPEPEVE ELESYWCEAA TATRSVRKRS PWLQLPGPGS PLDPASTTAP APWAAALAPG NRPLSFRKPP VSRSVPLVTS VRNTTSTGLQ FPASGAPTAG PPACAPPTPL EQSAGALKPD VDLLLREMQ LKGLLSRVVL ELKEPQALRE LRGTPTPTS HFAVSPGTPE TTPVES
Source	E.coli
Target Names	FCRLB
Protein Names	Recommended name: Fc receptor-like B Alternative name(s): Fc receptor homolog expressed in B-cells protein 2 Short name= FREB-2 Fc receptor-like and mucin-like protein 2 Fc receptor-like protein 2 Fc receptor-related protein Y
Expression Region	18-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.