



Recombinant Human PRAME family member 17 (PRAMEF17)

Product Code	CSB-EP018612HU-B
Abbreviation	PRAMEF17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VTA0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLQSPSRLL ELAQSSLRN QFLTIFILDE LPREVFPLMF MEASSMRHFE ALKLMVQAWP FLRLPLGSLM KTPHLETLQA VLKGLDTLLA QKLRPRRWKL QVLDLRDVDG NFWTIWVGAR ALSCSPEAMS KRQTVEDYPR TGEHQPLKVF IDLCQKESTL DECLSYLCRW IHYRRGLVHL CCNKVQNYSM PTSSFRNLLK RVYPDSIQEL EIKRKCSLNK TGKFAPYLSQ MSNLRKLFLA FGYDDELYVS GQQQFVPDLD CPFLCLYYPQ MLYIRKISNI KEHLEHLLRC LKNPLGTFIF CHAYLADQDM ECLSQYPSLS QLKELHLIHI LMWTTNLEPL GALLEKVAAT LEILTLKDCQ IQDSQLRVLL PALSRCSQLT TFYFRGNETS TNALKDLLCH TGGLSKLGLE LYPAPLECLD NRGHVNWEIL APIRAELMCT LREVRQPKRI FFGPIPCPSC GSWPSEKVDF HLCS
Source	E.coli
Target Names	PRAMEF17
Protein Names	Recommended name: PRAME family member 17
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.