



Recombinant Human PRAME family member 7 (PRAMEF7)

Product Code	CSB-MP018624HU
Abbreviation	PRAMEF7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VXH5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSIRAPPRL ELARQRLLRD QALAISTMEE LPRELFPTLF MEAFSRRRCE TLKTMVQAWP FTRLPLGSLM KSPHLESLKS VLEGVDVLLT QEVRPRQSKL QVLDLRNVDE NFCDIFSGAT ASFPEALSQK QTADNCPGTG RQQPFMVFID LCLKNRTLDE CLTHLLEWVK QKRGLLHVCC KELQVFGMPI HSIIEVLNMV ELDCIQEVEV CCPWELSTLV KFAPYLGQMR NLRKLVLFNI RASACIPPDN KGQFIARFTS QFLKLDYFQN LSMHSVSFLE GHLDQLLRCL QASLEMVVMV DCLLSESDLK HSWCPSIRQ LKELDLRGVT LTHFSPEPLT GLLEQVVATL QTLDEDCGI MDSQLSAILP VLSRCSQLST FSFCGNLISM AALENLLRHT VGLSKLSLEL YPAPLESYDT QGALCWGRFA ELGAELMKTP RDLRQPKIIV FCTVPCPRCG IRASYDLEPS HCLC
Source	Mammalian cell
Target Names	PRAMEF7
Protein Names	Recommended name: PRAME family member 7
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.