



# Recombinant Human PRAME family member 7 (PRAMEF7)

<b>Product Code</b>	CSB-YP018624HU
<b>Abbreviation</b>	PRAMEF7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5VXH5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSIRAPPRL ELARQRLLRD QALAISTMEE LPRELFPTLF MEAFSRRRCE TLKTMVQAWP FTRLPLGSLM KSPHLESLKS VLEGVDVLLT QEVRPRQSKL QVLDLRNVDE NFCDIFSGAT ASFPEALSQK QTADNCPGTG RQQPFMVFID LCLKNRTLDE CLTHLLEWVK QKRGLLHVCC KELQVFGMPI HSIIEVLNMV ELDCIQEVEV CCPWELSTLV KFAPYLGQMR NLRKLVLFNI RASACIPPDN KGQFIARFTS QFLKLDYFQN LSMHSVSFLE GHLDQLLRCL QASLEMVVMV DCLLSESDLK HLSWCPSIRQ LKELDLRGVT LTHFSPEPLT GLLEQVVATL QTLDEDCGI MDSQLSAILP VLSRCSQLST FSFCGNLISM AALENLLRHT VGLSKLSLEL YPAPLESYDT QGALCWGRFA ELGAELMKTP RDLRQPKIIV FCTVPCPRCG IRASYDLEPS HCLC
<b>Source</b>	Yeast
<b>Target Names</b>	PRAMEF7
<b>Protein Names</b>	Recommended name: PRAME family member 7
<b>Expression Region</b>	1-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.