



Recombinant Rat Mitogen-activated protein kinase 4 (Mapk4)

Product Code	CSB-YP013457RA
Abbreviation	Mapk4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63454
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	GCGGNGLVLS ATDSRACRKV AVKKIVLSDA RSMKHALREI KIIRRLDHDN IVKVYEVLPK KGSDLQGELF KFSVAYIVQE YMETDLACLL EQGTLTEEHA KLFMYQLLRG LKYIHSANVL HRDLKPANIF ISTDVLVLI GDFGLARIAD QHYSKGYLS EGLVTKWYRS PRLLLSPNNY TKAIDMWAAG CILAEMLTGK MLFAGAHELE QMQLLDTIP VVREEDKEEL LRVMPFSVSS TWEVKRPLRK LLPDVNREAI DFLEKILTFS PMDR
Source	Yeast
Target Names	Mapk4
Protein Names	Recommended name: Mitogen-activated protein kinase 4 Short name= MAP kinase 4 Short name= MAPK 4 EC= 2.7.11.24 Alternative name(s): Extracellular signal-regulated kinase 4 Short name= ERK-4 MAP kinase isoform p63
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Mitogen-activated protein kinase 4 is a member of the mitogen-activated protein kinase family. Tyrosine kinase growth factor receptors activate mitogen-activated protein kinases which then translocate into the nucleus where it phosphorylates nuclear targets.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.