



Recombinant Human ORC ubiquitin ligase 1 (OBI1)

Product Code	CSB-EP019878HU-B
Abbreviation	RNF219
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5W0B1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAQTVQNVTL SLTLPITCHI CLGKVRQPVI CINNHVFC SI CIDLWLKNNS QCPACRVPIT PENPCKEIIG GTSESEPMLS HTVRKHLRKT RLELLHKEYE DEIDCLQKEV EELKSKNLSL ESQIKTILDP LTLVQGNQNE DKHLVTDNPS KINPETVAEW KKKLRTANEI YEKVKDDVDK LKEANKKLLK ENGGLVRENL RLKAEVDNRS PQKFGRFAVA ALQSKVEQYE RETNRLKKAL ERSDKYIEEL ESQVAQLKNS SEEKEAMNSI CQTALSADGK GSKGSEEDVV SKNQGDSARK QPGSSTSSSS HLAKPSSSRL CDTSSARQES TSKADLNCSK NKDLYQE QVE VMLDVTDTSM DTYLEREWGN KPSDCVPYKD EELYDLPAPC TPLSLSLCLQL STPENRESSV VQAGGSKKHS NHLRKLVFDD FCDSSNVS NK DSSEDDISRS ENEKKSECFS SPKTGFWDCC STSYAQN LDF ESSEGNTIAN SVGEISSKLS EKSGLC LSKR LNSIRSFEMN RTRTSSEASM DAAYLDKISE LDSMMSESDN SKSPC NNGFK SLDLDGLSKS SQGSEFLEEP DKLEEKTELN LSKGSLTNDQ LENGSEWKPT SFFLLSPSDQ EMNEDFSLHS SSCPVTNEIK PPSCLFQTEF SQGILLSSSH RLFEDQRFGS SLFKMSSEM H SLHNHLQSPW STSFVPEKRN KNVNQSTKRK IQSSLSSASP SKATKS
Source	E.coli
Target Names	RNF219
Protein Names	Recommended name: RING finger protein 219
Expression Region	1-726
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.