



Recombinant Mouse Nucleolar protein 10 (No110)

Product Code	CSB-EP015918MO
Abbreviation	No110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RJG1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MQVSSLNEVK IYSLSCGKSL PEWLSDRKKR ALQKKNVDVR RRIELIQDFE MPTVCTTIKV SKDGQYILAT GTYKPRVRCY DTYQLSLKFE RCLDSEVVTF EILSDDYSKI VFLHNDRYIE FHSQSGFYK TRIPKFGDRDF SYHYPSCDLY FVGASSEVYR LNLEQGRYLN PLQTDAAENN VCDINAVHGL FATGTIEGRV ECWDPRVRKR VGVLDALNS VTADSEINSL PTISALKFNG ALSMAVGTST GQVLLYDLRS DKPLLVKDHQ YGLPIKSVHF QDSL DLVLSA DSRIVKMWNK DSGKIFTSLE PEHDLNDVCL YPSSGMLLTA NESPKMGIYY IPVLGPAPRW CSFLDNLTEE LEENPESTVY DDYKFVTKKD LENLGLTHLI GSPFLRAYMH GFFMDIRLYH KVKLMVNPFA YEEYRKDKIR QKIEETRAQR VQLKKLPKVN KELALKLIEE EEEKQKSTLK KVKSLPNIL TDDRFKVMFE NPDFQVDEES EEFRLNPLV SRISEKRKKQ LRLLEQQELK DEEEEEPEGK PSDAESSSESS DDEKGVVEEV RKQRLLQQE ERVKRQEQLK EDQQTVLKPKQ FYEIKAGEEF RSFKESATKQ RLMNKTLEDR LKLEAKHGTL NVSDTTVGSK QLTFTLKRSE QQKKQQEAEK LHRQERKNLR RSASHLRSRP RRGRPFH
Source	E.coli
Target Names	No110
Protein Names	Recommended name: Nucleolar protein 10
Expression Region	1-687
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.