



# Recombinant Mouse Nucleolar protein 10 (NoI10)

<b>Product Code</b>	CSB-YP015918MO
<b>Abbreviation</b>	NoI10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RJG1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQVSSLNEVK IYSLSCGKSL PEWLSDRKKR ALQKKNVDVR RRIELIQDFE MPTVCTTIKV SKDGQYILAT GTYKPRVRCY DTYQLSLKFE RCLDSEVVTF EILSDDYSKI VFLHNDRYIE FHSQSGFYK TRIPKFGGRDF SYHYPSCDLY FVGASSEVYR LNLEQGRYLN PLQTDAAENN VCDINAVHGL FATGTIEGRV ECWDPRVRKR VGVLDALNS VTADSEINSL PTISALKFNG ALSMAVGTST GQVLLYDLRS DKPLLVKDHQ YGLPIKSVHF QDSL DLVLSA DSRIVKMWNK DSGKIFTSLE PEHDLNDVCL YPSSGMLLTA NESPKMGIYY IPVLGPAPRW CSFLDNLTEE LEENPESTVY DDYKFVTKKD LENLGLTHLI GSPFLRAYMH GFFMDIRLYH KVKLMVNPFA YEEYRKDKIR QKIEETRAQR VQLKKLPKVN KELALKLIEE EEEKQKSTLK KVKSLPNIL TDDRFKVMFE NPDFQVDEES EEFLLNPLV SRISEKRKKQ LRLLEQQELK DEEEEEPEGK PSDAESSSESS DDEKGVVEEV RKQRLLQQE ERVKRQEQLK EDQQTVLKPKQ FYEIKAGEEF RSFKESATKQ RLMNKTLEDR LKLEAKHGTL NVSDTTVGSK QLTFTLKRSE QQKKQQEAEK LHRQERKNLR RSASHLRSRP RRGRPFH
<b>Source</b>	Yeast
<b>Target Names</b>	NoI10
<b>Protein Names</b>	Recommended name: Nucleolar protein 10
<b>Expression Region</b>	1-687
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.