



Recombinant Rat Proline-rich nuclear receptor coactivator 1 (Pnrc1)

Product Code	CSB-MP018289RA
Abbreviation	Pnrc1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63647
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MTVVSVPQRD LALGGRLAPL GSLGALPMVS PAPPPLPRLP DPRALPPTLF LPHFLGGDGP CLARASHRCT VQLQLGPGRC PTAAPKKRRK KKVRASPAGQ LPSRFHQFQQ HRPSLEGGRS PVPGPAAQE ERGPGATALL YRQPPLAKEV LKSKMGKSEK ISPPHSQLVH GIHLCEQPKI SRQKSKFNLP LTKITSAKRN ENDFWQDSAS SDRMQKQEKK PLKNTENIKN NHLKKSFLT EGSQKENYAG AKFSDPPSPS VLPKPPSHWM GSTTENPSQS RELMAVHLKT LLKVQT
Source	Mammalian cell
Target Names	Pnrc1
Protein Names	Recommended name: Proline-rich nuclear receptor coactivator 1 Alternative name(s): Proline-rich protein 2
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.