



Recombinant Human Zinc finger protein 626 (ZNF626)

Product Code	CSB-MP026904HU
Abbreviation	ZNF626
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68DY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGPLQFRDVA IEFSLLEWHC LDTAQRNLYR NVMLNYSNL VFLGITVSKP DLITCLEQGR KPLTMKRNEM IAKPSVMCSH FAQDLWPEQS MKDSFQKVVL RRYEKCEHDN LQLKKGCSV DECKVHKEGY NELNQCLTTT PRKICQCDKY VKVLHQFPNS NGQKRGHTGK KPFKYIECGK AFKQFSTLTT HKKIHTGGKP YKCEECGKAF NHSCSLTRHK KIHTGEKPYK CEECGKAFKH SSTLTTHKRN HTGEKPYKCD KCGKAFMSSS TSKHEIIHT EKKPYKCEEC GKAFNRSSTL TTHKIIHTGE KPYKCEECDK AFKYSYTLTT HKRIHTEDKP YKCEECGKAF KYSSTLTTHK RIHTGEKPYK CEECGKAFKR SSDLTTHKII HTGEKPYKCE ECGKAFKYSS NLTTHKKIHT GERPYKCEEC GKAFNQSSIL TTHRRIHTGE KFYKCEECGK AFKCSSNLTT HKKIHTGERP YKCEECGKAF NQSSILTHE RIILERNSTN VKNVAKPSSG PHTLLHIR
Source	Mammalian cell
Target Names	ZNF626
Protein Names	Recommended name: Zinc finger protein 626
Expression Region	1-528
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.