



Recombinant Human Dermokine (DMKN)

Product Code	CSB-BP006965HU
Abbreviation	DMKN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6E0U4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GPLQSGEES TGTNIGEALG HGLGDALSEG VGKAIGKEAG GAAGSKVSEA LGQGTREAVG TGVRQVPFGF VADALGNRVG EAAHALGNTG HEIGRQAEDV IRHGADAVRG SWQGVPGHNG AWETSGGHGI FGSQGGLGGQ GQGNPGLGT PWVHGYPGNS AGSFGMNPQG APWGQGGNGG PPNFGTNTQG AVAQPGYGSV RASNQNEGCT NPPPSGSGGG SSNSGGGSGS QSGSSGSGSN GDNNNGSSSG GSSSGSSSGG SSGGSSGGSS GNSGGSRGDS GSESSWGSST GSSSGNHGGS GGGNGHKPGC EKPGNEARGS GESGIQNSET SPGMFNFDTF WKNFKSKLGF INWDAINKNQ VPPPSTRALL YFSRLWEDFK QNTPLFNWKA IIEGADASSL QKRAGRDDQN YNYNQHAYPT AYGGKYSVKT PAKGGVSPSS SASRVQPGLL QWVKFW
Source	Baculovirus
Target Names	DMKN
Protein Names	Recommended name: Dermokine Alternative name(s): Epidermis-specific secreted protein SK30/SK89
Expression Region	22-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene is upregulated in inflammatory diseases, and it was first observed as expressed in the differentiated layers of skin. The most interesting aspect of this gene is the differential use of promoters and terminators to generate isoforms with unique cellular distributions and domain compliments. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene, but the full-length nature of some of them has not been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.