



# Recombinant Pongo abelii PHD finger protein 6 (PHF6)

<b>Product Code</b>	CSB-EP017917PYX
<b>Abbreviation</b>	PHF6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R5Z2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSSVEQKKG PTRQRKCGFC KSNRDKECGQ LLISENQKVA AHHKCMLFSS ALVSSHSDNE SLGGFSIEDV QKEIKRGTKL MCSLCHCPGA TIGCDVKTCH RTYHYHCALH DKAQIREKPS QGIYMYCRK HKKTAHNSEA DLEESFNEHE LEPSSPKSKK KSRKGRPRKT NFKGLSEDTR STSSHGTDEM ESSSYRDRSP HRSSPSDTRP KCGFCHVGEE ENEARGKLHI FNAKAAAHY KCMLFSSGTV QLTTTSRAEF GDFDIKTVLQ EIKRGKRMKC TLCSQPGATI GCEIKACVKT YHYHCGVQDK AKYIENMSRG IYKLYCKNHS GNDRDEEDE ERESKSRGKV EIDQQQLTQQ QLNGN
<b>Source</b>	E.coli
<b>Target Names</b>	PHF6
<b>Protein Names</b>	Recommended name: PHD finger protein 6
<b>Expression Region</b>	2-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.