



Recombinant Human Zinc finger protein 684 (ZNF684)

Product Code	CSB-MP026950HU
Abbreviation	ZNF684
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T5D7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MISFQESVTF QDVAVDFTAE EWQLLDCAER TLYWDVMLEN YRNLISVGCP ITKTKVILKV EQGQEPWMVE GANPHESSPE SDYPLVDEPG KHRESKDNFL KSVLLTFNKI LTMERIHYN MSTSLNPMRK KSYKSFEKCL PPNLDLLKYN RSYTVENAYE CSECGKAFKK KFHFIRHEKN HTRKKPFECN DCGKAYSRA HLATHQKIHN GERPFVCNDC GKAFMHKAQL VVHQRLLHTGE KPYECSQCGK TFTWNSSFNQ HVKSHTLEKS FECKEKGKTF RYSSSLYKHS RFHTGEKPYQ CIICGKAFGN TSVLVTHQRI HTGEKPYSCI ECGKAFIKKS HLLRHQITHT GEKPYECNRC GKAFSQKSNL IVHQKIHT
Source	Mammalian cell
Target Names	ZNF684
Protein Names	Recommended name: Zinc finger protein 684
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.