



Recombinant Mesocricetus auratus Synaptonemal complex protein 3 (SYCP3)

Product Code	CSB-EP022993MRG-B
Abbreviation	SYCP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60547
Product Type	Recombinant Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Purity	>85% (SDS-PAGE)
Sequence	MVPGGRKHSG KSGKPPLVDQ AKTAFDFEKE DKELSGSEED VADEKTPVID KHGKKRSGGL VEDVGGEVQN MLEKFGADIN KALLAKRKRI EMYTKASFKA SNQKIEQIWK TQQEEIQKLN SEYSQQFMSV LQQWELDMQK FEEQGEKLTN LFRQQQKIFQ QSRIVQSQR L KAIKQLHEQF IKNLEDVEKN NDNLFTGTQS ELKKEMAMLQ KKVMMETQQQ EMANVRKSLQ SMLF
Source	E.coli
Target Names	SYCP3
Protein Names	Recommended name: Synaptonemal complex protein 3 Short name= SCP-3 Alternative name(s): Meiotic chromosome core protein Synaptonemal complex protein COR1
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.