




Recombinant Mouse Mesothelin (Msln), partial

Product Code	CSB-EP015044MO
Relevance	Mbrane-anchored forms may play a role in cellular adhesion.
Abbreviation	Recombinant Mouse Msln protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61468
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	DAEQKACPPGKEPYKVDEDLIFYQNWELEACVDGTMLARQMDLVNEIPFTYE QLSIFKHKLDKTYPQGYPELIQQLGHFFRYVSPEDIHQWNVTSPTVKTLLKV SKGQKMNAQAIALVACYLRGGGQLDEDMVKALGDIPLSYLCDFSPQDLHSVP SSVMWLVGPDLDKCSQRHLGLLYQKACSAFQNVSGLEYFEKIKTFLGGASV KDLRALSQHNVSMDIATFKRLQVDSLVLGLSVAEVQKLLGPNIVDLKTEEDKSPV RDWLFQRHQKDLDRGLGLQGGIPNGYLVLDNFVREAFS
Source	E.coli
Target Names	Msln
Expression Region	298-600aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	50.1kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.