



Recombinant Danio rerio Nuclear prelamin A recognition factor (narf)

Product Code	CSB-BP015454DIL
Abbreviation	narf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHP6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSEVKIGRKK EKCENCTKQC NKKQSDESIN SLQEKGVASG EVRESTSVQQ VLLSACLSCD GCISEDEGKR ISQQNLDEIN HVLALNKKCD TSKHKILVVS VCPQSVFFA VKFQLDVSA AQLKCGFLKS VGVHYVFDTT IAASFSILES QREFVQRYRR KHHDANAMPM FTSSCPGWIR YAERVLGSVV TPHICTARSP QQIMGSLVKN FFARQQKLTPEQVFHVVVAP CFDKKLEAVR DEFYNSILES RDVDCVLTSG EILHLMEQSK VTVEEVDSAP LDHVICEISD SGLTRHDGRG SEGFLEHIFK HAAKEIFGLD VQDIVYKTLR NRDFQEVALE RDGETLLQFA AVYGFRNIQT LVHRMRKGKV PYQLIEVLSC PGGCLSGRGQ AEGEGGRPDR TLVQQMEEHY SSLPVRLPET NPEVQRLYQD WLDGHDSPPHA QQSLHTQYKD HTHLPAHITN PDIQW
Source	Baculovirus
Target Names	narf
Protein Names	Recommended name: Nuclear prelamin A recognition factor
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.