



Recombinant Human Ventral anterior homeobox 1 (VAX1)

Product Code	CSB-MP025806HU
Abbreviation	VAX1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SQQ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFGKPKMDV RCHSDAEAAR VSKNAHKESR ESKGAEGNLP AAFLKEPQGA FSASGAAEDC NKSNSNSAAD PDYCRRILVR DAKGSIREII LPKGLDLDRP KRTRTSFTAE QLYRLEMETFQ RCQYVVGRER TELARQLNLS ETQVKVWFQN RRTKQKKDQG KDELRSVVS ETAATCSVLR LLEQGRLLS PGLPALLPPC ATGALGSALR GPSPALGAG AAAGSAAAAA AAPGPAGAA SPHPPAVGGA PGPAGPAGGG LHAGAPAAGH SLFSLPVPVSL LGSVASRLSS APLTMAGSLA GNLQELSARY LSSSAFEPYS RTNNKEGAEK KALD
Source	Mammalian cell
Target Names	VAX1
Protein Names	Recommended name: Ventral anterior homeobox 1
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a homeo-domain containing protein from a class of homeobox transcription factors which are conserved in vertebrates. Genes of this family are involved in the regulation of body development and morphogenesis. The most conserved genes, called HOX genes are found in special gene clusters. This gene belongs to the VAX subfamily and lies in the vicinity of the EMX homeobox gene family. Another member of VAX family is located on chromosome 2. The encoded protein may play an important role in the development of anterior ventral forebrain and visual system. Multiple transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.