



Recombinant Human PHD finger protein 19 (PHF19)

Product Code	CSB-EP017907HU-B
Abbreviation	PHF19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T6S3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MENRALDPGT RDSYGATSHL PNKGALAKVK NNFKDLMSKL TEGQYVLCRW TDGLYYLGKI KRVSSSKQSC LVTFEDNSKY WVLWKDIQHA GVPGEETPCN ICLGKTSGPL NEILICGKCG LGYHQQCHIP IAGSADQPLL TPWFRCRCIF ALAVRKGGAL KKGAIARTLQ AVKMWLSYQP EELEWDSPHR TNQQQCYYCYC GPGGEWYLRM LQCYRCRQWF HEACTQCLNE PMMFGDRFYL FFCVCNQQGP EYIERLPLRW VDVVHLALYN LGVQSKKKYF DFEEILAFVN HHWELLQLGK LTSTPVTDRG PHLLNALNSY KSRFLCGKEI KKKKCIFRLR IRVPPNPPGK LLPDKGLLPN ENSASSELRK RGKSKPGLLP HEFQQQKRRV YRRKRSKFLLEDAIPSSDFT SAWSTNHHLA SIFDFTLDEI QSLKSASSGQ TFFSDVDSTD AASTSGSAST SLSYDSRWTV GSRKRKLA AKAYMPLRAKRW AAELDGRCPSSDSSAEGASVP ERPDEGIDSH TFESISSEDDSSLSHLKSSIT NYFGAAGRLA CGEKYQVLAR RVTPEGKVQY LVEWEGTTPY
Source	E.coli
Target Names	PHF19
Protein Names	Recommended name: PHD finger protein 19 Alternative name(s): Polycomb-like protein 3 Short name= hPCL3
Expression Region	1-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.