



Recombinant Danio rerio Centromere protein L (cenpl)

Product Code	CSB-YP005215DIL
Abbreviation	cenpl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DRD4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MFNLLPAFSL GVLLIVMEGR GSAVNTPASR PVTRKSKSYG RSCVEAPSYF WQTPGHLTVL KVPTSRGAAK SQIITNKVDP EHIALLVKNE WKLSYVTPH HFRHTQLKSY AKHLSAFIVA EKQQGVAVEV GLDAGFKVTF TAVLGLAETD EDAETVFIQI QSKPVFAAAA DAPKVVWRGW LTCVNGDPEY LRSLPPDFVS LPLFCTSGPE SLTNLVKSWF ERAFDCNFGS LALNSSTLNW LAALWTGCHP ACNIRYLKLS WSLPTEPPLD VLYTVNPQDA WELWNSIHPQ ENTDDRDIK DVQAFMSGME THFFRHFKIH LSAGTLLKVS TALGSAHHDG KIKIGNSDYI TTLLALLTEC ALLKMPT
Source	Yeast
Target Names	cenpl
Protein Names	Recommended name: Centromere protein L Short name= CENP-L
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.