



Recombinant Human Secreted frizzled-related protein 4 (SFRP4)

Product Code	CSB-MP021140HU
Abbreviation	SFRP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6FHJ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	VR GAPCEAVRIP MCRHMPWNIT RMPNHLHHST QENAILAIEQ YEELVDVNCS AVLRFFLCAM YAPICTLEFL HDPIKPCKSV CQRARDDCEP LMKMYNHSWP ESLACDELPV YDRGVCISPE AIVTDLPELV KWIDITPDMM VQERPLDVDC KRLSPDRCKC KVKKPTLATY LSKNYSYVIH AKIKAVQRSG CNEVTTVVDV KEIFKSSSPI PRTQVPLITN SSCQCPHILP HQDVLIMCYE WRSRMMLLEN CLVEKWRDQL SKRSIQWEER LQEQRRTVQD KKKTAGRISR SNPPKPKGKP PAKPASPCK NIKTRSAQKR TNPKRK
Source	Mammalian cell
Target Names	SFRP4
Protein Names	Recommended name: Secreted frizzled-related protein 4 Short name= sFRP-4 Alternative name(s): Frizzled protein, human endometrium Short name= FrpHE
Expression Region	19-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.