



Recombinant Human INO80 complex subunit C (INO80C)

Product Code	CSB-MP011726HU
Abbreviation	INO80C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PI98
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAQIPIVAT TSTPGIVRNS KKRPAASPHN GSSGGGYGAS KKKKASASSF AQGISM EAMS ENKMVPSEFS TGPVEKAAKP LPFKDPNFVH SGHGGAVAGK KNRTWKNLKQ ILASERALPW QLNDPNYFSI DAPPSFKPAK KYSDVSGLLA NYTDPQSKLR FSTIEEFSYI RRLPSDVVTG YLALRKATSI VP
Source	Mammalian cell
Target Names	INO80C
Protein Names	Recommended name: INO80 complex subunit C Alternative name(s): IES6 homolog Short name= hles6
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.