



Recombinant Human Docking protein 3 (DOK3)

Product Code	CSB-YP007109HU
Abbreviation	DOK3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L591
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTRGARLRSD ARAQLNQLSL DGGTGSGQKG KCEEFSSLS SVSPGLEAAA LLAVTMDPL ETPIKDGILY QQHVKFGKCC WRKVVALLYA GGPSGVARLE SWEVRDGGGLG AAGDRSAGPG RRGERRVIRL ADCVSVLPAD GESCPRTDGA FLLTTTTERSH LLAAQHRQAW MGPICQLAFP GTGEASSGST DAQSPKRGLV PMEENSIYSS WQEVGEFPVV VQRTEAATRC QLKGPALLVL GPDAILQREA KGTQALYSWP YHFLRKFGSD KGVFSFEAGR RCHSGEGLFA FSTPCAPDLC RAVAGAIARQ RERLPELTRP QPCPLPRATS LPSLDTPGEL REMPGPPEPP TSRKMHLAEP GPQSLPLLLG PEPNDLASGL YASVCKRASG PPGNEHLYEN LCVLEASPTL HGGEPEPHEG PGSRSPTTSP IYHNGQDLSW PGPANDSTLE AQYRRLLELD QVEGTGRPDP QAGFKAKLVT LLSRERRKGP APCDRP
Source	Yeast
Target Names	DOK3
Protein Names	Recommended name: Docking protein 3 Alternative name(s): Downstream of tyrosine kinase 3
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.