



Recombinant Human Zinc finger protein 530 (ZNF530)

Product Code	CSB-BP026818HU
Abbreviation	ZNF530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P9A1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAALRAPDQ QVFVAFEDVA IYFSQEEWEL LDEMQRLLYR DVMLLENFAVM ASLGCWCGAV DEGTPSAESV SVEELSQGRT PKADTSTDKS HPCEICTPVL RDILQMIELH ASPCGQKLYL GGASRDFWMS SNLHQLQKLD NGEKLFKVDG DQASFMMNCR FHVSGKPFTF GEVGRDFSAT SGLLQHQVTP TIERPHSRIR HLRVPTGRKP LKYTESRKSF REKSVFIQHQ RADSGERPYPK CSECGKSFSQ SSGFLRHRKA HGRTRTHECS ECGKSFSRKT HLTQHQRVHT GERPYDCSEC GKSFRQVSVL IQHQRVHTGE RPYECSECGK SFSHSTNLYR HRSHTSTRP YECSECGKSF SHSTNLFRHW RVHTGVRPVE CSECGKAFSC NIYLIHHQRF HTGERPYVCS ECGKSGGQKS VLIQHQRVHT GERPYECSEC GKVFSQSSGL FRHRAHTKT KPYECSECEK SFCKTDLIR HQTVHTGERP YECVCGKSF IRKTHLIRHQ TVHTNERPVE CDECGKSYSQ SSALLQHRRV HTGERPYECP ECGKSFTRKN HLIQHKTVHT GERPYECSEC GKSFSQSSGL LRHRRVHVQ
Source	Baculovirus
Target Names	ZNF530
Protein Names	Recommended name: Zinc finger protein 530
Expression Region	1-599
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.