



# Recombinant Human RING finger protein 207 (RNF207)

<b>Product Code</b>	CSB-BP019870HU
<b>Abbreviation</b>	RNF207
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZRF8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSGAIFGPLE GPSSLDAPSI HPLVCPLCHV QYERPCLLDC FHDFCAGCLR GRATDGRLTC PLCQHQTVLK GPSGLPPVDR LLQFLVDSSG DGVEAVRCAN CDLECSEQDV ETTYFCNTCG QPLCARCRDE THRARMFARH DIVALGQRSR DVPQKCTLHA EPYLLFSTDK KLLLCIRCFR DMQKESRAHC VDLESAYVQG CERLEQAVLA VKALQTATRE AIALQAMVE EVRHSAEEEE DAIHALFGSM QDRLAERKAL LLQAVQSQYE EKDKAFKEQL SHLATLLPTL QVHLVICSSF LSLANKAEFL DLGYELMERL QGIVTRPHHL RPIQSSKIAS DHRAEFARCL EPLLLLGP RR VAAAASGANT LAGGLGPKAL TGPHCPSPVG KMSGSPVQKP TLHRSISTKV LLAEGENTPF AEHCRHYEDS YRHLQAEMQS LKDQVQELHR DLTKHHS LIK AEIMGDVLHK SLQLDVQIAS EHASLEGMRV VFQEIWEEAY QRVANEQEIY EAQLHDLLQL RQENAYLTTI TKQITPYVRS IAKVKERLEP RFQAPVDEQS ESLQNTHDDS RNNAASARNN PGSVPEKREK TSEPKGNSWA PNGLSEEPLL KNMDHHR SKQ KNGGDVPTWR EHPT
<b>Source</b>	Baculovirus
<b>Target Names</b>	RNF207
<b>Protein Names</b>	Recommended name: RING finger protein 207
<b>Expression Region</b>	1-634
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.