



Recombinant Human Speckle-type POZ protein-like (SPOPL)

Product Code	CSB-BP022602HU
Abbreviation	SPOPL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IQ16
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSREPTPLP GDMSTGPIAE SWCYTQVKVV KFSYMWTTINN FSFCREEMGE VLKSSTFSSG PSDKMKWCLR VNPKGLDDES KDYLSTLYLLL VSCPKEVRA KFKFLLNAK REETKAMESQ RAYRFVQGKD WGFKKFIRRD FLLDEANGLL PDDKLTLFCE VSVVQDSVNI SGHTNTNTLK VPECRLAEDL GNLWENTRFT DCSFFVRGQE FKAHKSVLAA RSPVFNAMFE HEMEESKKNR VEINDLDPEV FKEMMRFIYT GRAPNLDKMA DNLLAAADKY ALERLKVMECE EALCSNLSVE NVADTLVLAD LHSAEQLKAQ AIDFINRCSV LRQLGCKDGK NWNSNQATDI METSGWKSMI QSHPHLVAEA FRALASAQCP QFGIPRKRLK QS
Source	Baculovirus
Target Names	SPOPL
Protein Names	Recommended name: Speckle-type POZ protein-like Alternative name(s): HIB homolog 2 Roadkill homolog 2
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.