



Recombinant Human Zinc finger protein 774 (ZNF774)

Product Code	CSB-BP027030HU
Abbreviation	ZNF774
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NX45
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MWLGTSGKSG LPGHCLENPL QECHPAQLEE WALKGISRPS VISQPEQKEE PWVLPQNF E ARKIPRESHT DCEHQVAKLN QDNSETAEQC GTSSERTNKD LSHTLSWGGN WEQGLELEGQ HGTLPGEGQL ESFSQERDLN KLLDGYVGEK PMCAECGKSF NQSSYLIRHL RTHTGERPYT CIECGKGFQK SSDLVTHRRRT HTGEKPYQCK GCEKKFSDSS TLIKHQRTHT GERPYECPEC GKTFGRKPHL IMHQRTHTGE KPYACLECHK SFSRSSNFIT HQRTHTGVPK YRCNDCGESF SQSSDLIKHQ RTHTGERPFK CPECCKGFRD SSHFVAHMST HSGERPFSCP DCHKSFSQSS HLVTHQRTHT GERPFKCENC GKGADSSAL IKHQRIHTGE RPYKCGECGK SFNQSSHFIT HQRHILGDRP YRCPECCKTF NQRSHFLTHQ RTHTGEKPFH CSKCNKSFRQ KAHLLCHQNT HLI
Source	Baculovirus
Target Names	ZNF774
Protein Names	Recommended name: Zinc finger protein 774
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.