



Recombinant Human Carboxylesterase 3 (CES3)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP005260HU |
| Abbreviation | CES3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6UWW8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | TGPE VAQPEVDTTL GRVRGRQVGV KGTDRLVNVF LGIPFAQPPL GPDRFSAPHP AQPWEGVRDA STAPPMCLQD VESMNSSRFV LNGKQQIFSV SEDCLVLNVY SPAEVPAGSG RPYMVVWHGG ALITGAATSY DGSALAAAYGD VVVTVQYRL GVLGFFSTGD EHAPGNQGFL DVVAALRWVQ ENIAPFGGDL NCVTVFGGSA GGSIIISGLVL SPVAAGLFHR AITQSGVITT PGIIDSHWP LAQKIANLA C SSSSPAEMV QCLQQKEGEE LVLSKCLKNT IYPLTVDGTV FPKSPKELLK EKP FHSVPFL MGVNNHEFSW LIPRGWLLD TMEQMSREDM LAISTPVLTS LDVPPMMPT VIDEYLG SNS DAQAKCQAFQ EFMGDVFINV PTVSFSRYLR DSGSPVFFYE FQHRPSSFAK IKPAWVKADH GAEGAFVFGG PFLMDESSL AFPEATEEEK QLSLTMMAQW THFARTGDPN SKALPPWPQF NQAQYLEIN PVPRAGQKFR EAWMQFWSET LPSKIQQWHQ KQKNRKAQED L |
| Source | Mammalian cell |
| Target Names | CES3 |
| Protein Names | Recommended name: Carboxylesterase 3 EC= 3.1.1.1 Alternative name(s): Liver carboxylesterase 31 homolog |
| Expression Region | 27-571 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | Carboxylesterase 3 is a member of a large multigene family. The enzymes encoded by these genes are responsible for the hydrolysis of ester- and amide-bond-containing drugs such as cocaine and heroin. They also hydrolyze long-chain fatty acid esters and thioesters. The specific function of this enzyme has not yet been determined; however, it is speculated that carboxylesterases may play a role in lipid metabolism and/or the blood-brain barrier system. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.